

## • Concentric Tubing Masks

Using Concentric Tubing design, sturdy and durable, can be placed arbitrarily on the desktop;  
The entire mask (including tubing) is made of non-metallic material suitable for MRI environments;

- The mask can be fixed by fixed bracket (68625) or the fixed platform (68626), height of the mask is adjustable, and it can also be fixed and placed in MRI, PET and SPECT imaging systems



- Simultaneous experiments can be performed on 2-5 rats or mice with five channels of anesthesia mask



- Universal bracket (R510-55) can be used to fix the R510-32, R510-33 Tubing mask.



### Order Information:

Cat No.	Product Description	Cat No.	Product Description
68625	MRI Compatible Mounting Bracket for Tubing or Cone Mask, optional	68626	MRI Compatible Mounting Platform for Tubing or Cone Mask-rat/mouse
R510-55	Universal Stand for Tube Anesthesia Mask		