

Induction Chambers

For rapid induction of anesthesia before animal surgery, takes only 2-5 minutes; Animal status can be observed in real time during anesthesia.

The anesthetic gas is inducted into the chamber, and exhausted to the connected charcoal canister. With following features:

- Transparent plexiglass material, durable, easy to clean, convenient to observe the anesthetic state of animals;
- Positive pressure gaskets at lids;
- Both ports on one side to reduce tubing in procedure area.



Order Information:

Cat No.	Product Description	Remarks
V100	Induction Chamber- mouse (15cm*10cm*10cm)	
V101	Induction Chamber- rat (24cm*12cm*18cm)	Contains 1 plexiglass box, 1.2m bellows and 1 gas filter canister (R510-31S)
V102	Induction Chamber- rabbit/cat (40cm*18.5cm*25cm)	
V103	Induction Chamber- large animal (50cm*30cm*30cm)	