

Injection & Blood Collection Needles

For Laboratory Animals size 25G with special port access






Can be Used For Both Injection and Blood Collection






- Connects to hematocrit tubes or syringes and are adaptable for use in various treatments.
- Connects to both plastic and glass hematocrit tubes of various widths.

Enters Blood Vessels Easily





- Wings allow stable grip.
- Inserts approximately parallel to the animals' tail without being angled.
- Higher rate of successful venipuncture that lead to improvement of experiment precision.

Minimizes Damage

- as needles are made to be non-invasive on blood vessels and on the body, injection and blood collection can be performed frequently.
- The closed system protects the safety of animals by enabling direct blood collection from the caudal vein.
- Bood is collected through the natural flow of capillaries which keeps damage to red blood cells minimal and prevents hemolysis.

