



Aniphy

Multi-output flowmeter



R510-07



R510-09

R510-07 five-output flowmeter and R510-09 can be used in conjunction with R520 and R540 series small animal anesthesia machines to meet the simultaneous anesthesia of 5 or 2 animals, and each port can adjust the flow rate independently.

Extended flowmeter



RWD anesthesia is equipped with an extended flow meter mounting position, which expands a variety of gas flowmeters to meet the experimental needs of a variety of gas supply.

Order Information:

Cat No.	Product Description	Remarks
R510-07	Five-output flowmeter - range 1L / min	Standard contains 5.3 meter silicon tube (ID:6.0mm)
R510-09	Dual Channel Flowmeter - Range 0.1-1L / min	Standard contains 2.2 meter silicon tube (ID:6.0mm)
R500-N2O-1E	Expansion Assembly with 0.1-1LPM N ₂ O Flowmeter	Silicone tube and Y-three way and other accessories
R600-O2-4E	Expansion Assembly with 0.1-4LPM O ₂ Flowmeter	Silicone tube and Y-three way and other accessories