


## ● Cone Masks

- Special designed with inlet tubing inside make the animal anesthesia exhaust gas can be absorbed by gas filter canister;
- Allows access to eyes, ears and top of head;
- The cone masks can be installed onto 68630 operation bed for single cone mask (300\*210\*75mm) or 68657 manifold for five cone masks



### Order Information:

Cat No.	Product Description	OD (size)	ID(size)	D(size)
68653	Cone Mask with Tubing for Mouse or Rat (<70g)-Tube	13mm	5mm	24mm
68654	Cone Mask with Tubing for Rat (150-250g)-Tube	20mm	10mm	28.5mm
68655	Cone Mask with Tubing for Rat (250-350g)-Tube	25mm	15mm	38mm
R-CT-6/15	Concentric Tubing, 6mm Inlet/15mm Outlet, 45cm			
68656	Manifold for Single Cone Mask			
68657	Manifold for Five Cone Masks			
68630	Anesthesia Operation Platform			
68658	Cone Mask for Mouse or Rat (<70g)	13mm	5mm	24mm
68659	Cone Mask for Rat(150-250g)	20mm	10mm	28.5mm
68660	Cone Mask for Rat(250-350g)	25mm	15mm	38mm