




## Small Animal Anesthesia Machines ▶

### R500 Compact Anesthesia Machine

Regular assembly (R500-A), vaporizer(R580), Inlet/Outlet End Caps(R580-IO) and tube accessory kit are included in the R500 anesthesia machine. Different air supply, induction chambers, small animal masks, filter canisters and gas evacuation apparatuses can be selected for anesthesia recovery of different animals.

**Suitable for rats, mice, hamsters, guinea pigs, rabbits, cats and other animals within 7kg**



Note: Oxygen flush is commonly used to flush out anesthetic gas to the filter canister before opening the induction chamber to avoid the escape of anesthetic gas into the personnel.

#### Order Information:

Cat No.	Product Description	Configuration	Remarks
R500IE	Compact Small Animal Anesthesia Device-Isoflurane, Easy Fill, 0-1L/0-4L	●	
R500IP	Compact Small Animal Anesthesia Device-Isoflurane, Pour Fill, 0-1L/0-4L	●	Contains connection tubing for mask&induction chamber
R500SE	Compact Small Animal Anesthesia Device-Sevoflurane, Easy Fill, 0-1L/0-4L	●	
R500SP	Compact Small Animal Anesthesia Device-Sevoflurane, Pour Fill, 0-1L/0-4L	●	
R510-S	Portable Stand for R510 Series Anesthesia Machine	○	Upgrade to portable type
R600-N2O-4E	Expansion Assembly with 0.1-4LPM N <sub>2</sub> O Flowmeter	○	
R510-09	2ch Flowmeter for Anesthesia, range 0.1-1.0L/min each	○	Simultaneous anesthesia for two animals ( adjusted independently )

Note: ●Recommended configuration ○ Optional