


Mice and Rats restraint devices



Greatly Contributes to
“Refinement”
of Animal Use
Alternatives



Bell Shaped Entrance

Designed for
smooth entry that
prevents injuries
and stress

- Minimizes animals' distress, and significantly improves efficiency
- Ideal for procedures on the tail
- Safe to sterilize in autoclaves

Product	Animal	Dimension (mm) Inner
Restraint ϕ 28	Mice	Diameter 28 \times L110,
Restraint ϕ 32	Mice	Diameter 32 \times L110,
Restraint ϕ 35	Mice	Diameter 35 \times L126,
Restraint ϕ 47	Rats	Diameter 47 \times L176,
Restraint ϕ 64	Rats	Diameter 64 \times L215,
Restraint ϕ 70	Rats	Diameter 70 \times L215,

 **Aniphy**
Instruments That Work !

INSTRUMENTS AND CONSUMABLE FOR THE
RESEARCH ON NEUROSCIENCES,
PHARMACOLOGY AND PHYSIOLOGY IN VIVO

Aniphy.com

+33 484 497 005

contact@aniphy.com